

Certificate Of Analysis

Product KAPA LTP Lib Prep Kit (8rxn)
 Kit Code KK8230
 Material Number 07961863001
 Lot Number 004859-1-1
 Expiry Date 2018-12-22

Components		Lot Number
KB8200	10X End-repair Buffer	00067368
KB8210	10X A-Tailing Buffer	00066708
KB8220	5X Ligation Buffer	00067732
KB8232	PEG/NaCl	00068466
KE8200	End-repair Enzyme Mix	00068754
KE8210	A-Tailing Enzyme	00067386
KE8220	DNA Ligase	00066922
KM2604	2X HiFi HotStart ReadyMix	00067494
KP8202	Illumina Lib Amp Primers	00068592

Parameter	Result
Contamination Assays	Conforms
Bacterial	Pass
Endonucleases	Pass
Exonucleases	Pass
Human	Pass
Library fragment	Pass
Functional Test	Conforms
Functional	Pass

This Lot meets the specification.

CofA auto-generated on 2016/12/28