

Certificate Of Analysis

Product KAPA HyperPlus (8rxn)
 Kit Code KK8510
 Material Number 07962380001
 Lot Number 004929-11-1
 Expiry Date 2017-11-24

Components		Lot Number
KB8500	End Repair & A-Tailing Buffer	00069327
KB8501	Ligation Buffer	00069331
KE8500	End-Repair and A-Tailing Enzyme	00069339
KE8501	DNA Ligase Enzyme	00069347
KM2604	2X HiFi HotStart ReadyMix	00067494
KP8202	Illumina Lib Amp Primers	00068592
KB8601	10X KAPA Frag Buffer	00069145
KE8601	5X KAPA Frag Enzyme	00068019
KB8600	Frag Conditioning Solution	00068741

Parameter	Result
Contamination Assays	Conforms
Bacterial	Pass
Endonucleases	Pass
Exonucleases	Pass
Human	Pass
Library fragment	Pass
Functional Test	Conforms
Functional	Pass

This Lot meets the specification.

CofA auto-generated on 2017/03/10