

Certificate Of Analysis

Product KAPA HyperPlus (96rxn)
 Kit Code KK8514
 Material Number 07962428001
 Lot Number 004654-3-1
 Expiry Date 2017-09-10

Components		Lot Number
KB8504	End Repair & A-Tailing Buffer	00067626
KB8505	Ligation Buffer	00067622
KE8504	End-Repair and A-Tailing Enzyme	00067606
KE8505	DNA Ligase Enzyme	00067870
KM2612	2X HiFi HotStart ReadyMix	00067534
KP8201	Illumina Lib Amp Primers	00067898
KB8603	10X KAPA Frag Buffer	00067996
KE8603	5X KAPA Frag Enzyme	00068012
KB8600	Frag Conditioning Solution	00067896

Parameter	Result
<u>Contamination Assays</u>	Conforms
Bacterial	Pass
Endonucleases	Pass
Exonucleases	Pass
Human	Pass
Library fragment	Pass
<u>Functional Test</u>	Conforms
Functional	Pass

This Lot meets the specification.

CofA auto-generated on 2016/11/11