

Certificate Of Analysis

Product KAPA LTP Lib Prep Kit PCR-free (48rxn)
Kit Code KK8233
Material Number 07961898001
Lot Number 005010-10-1
Expiry Date 2019-04-18

Components		Lot Number
KB8201	10X End-repair Buffer	00068729
KB8211	10X A-Tailing Buffer	00068725
KB8221	5X Ligation Buffer	00068710
KB8231	PEG/NaCl	00068304
KE8201	End-repair Enzyme Mix	00068755
KE8211	A-Tailing Enzyme	00068745
KE8221	DNA Ligase	00067872

Parameter		Result
Contamination Assays		Conforms
	Bacterial	Pass
	Endonucleases	Pass
	Exonucleases	Pass
	Human	Pass
	Library fragment	Pass
Functional Test		Conforms
	Functional	Pass

This Lot meets the specification.

CofA auto-generated on 2017/02/22